

## Volumetric KF titrators – tabular



Technical specification	870 KF Titrino Plus	890 Titrando	841 Titrando
Volumetric Karl Fischer titration (KFT)	yes	yes	yes
Titration to preset pH/mV endpoint (SET)	–	–	yes
pH calibration and pH measurement (CAL, MEAS)	–	–	yes
Real-time curve display	yes	yes	yes
Predefined methods	yes	yes	yes
Parallel titration	–	–	yes ( <i>tiamo</i> )
<b>Inputs, interfaces</b>			
Second, galvanically separated measuring interface	–	–	option
Automated reagent exchange with Dosino	–	yes	yes
MSB-port for connection of further Dosinos (addition of auxiliary solution,...)	–	3	4
MSB-port for connection of further Dosinos (Titration)	–	–	4
Attachment of stirrer	yes	yes	yes
Sample changer connection	–	yes	yes
Attachment of printer, balance, ...	USB (1)	USB (2)	USB (2)
Optional interfaces	RS 232C	RS 232C	RS 232C
<b>Software</b>			
<i>tiamo</i> <sup>™</sup>	–	yes	yes
PC Control or Touch Control	–	yes	yes